

**SafetyMate™ is the first smart device to verbally coach users through a medical emergency until professional help arrives.**



SafetyMate is the leading Emergency Instruction Device (EID), providing quick and efficient coaching during a medical emergency until professional help arrives. The portable EID provides immediate, easy-to-follow verbal instructions for emergency medical care. The user interacts with the device by responding to Yes/No questions, and SafetyMate uses the answers to provide step-by-step instructions for how to provide immediate and effective care. The protocols were developed in close cooperation with the American Red Cross and are based on the same information used by 911 operators.



**SAFETY + MATE™**

WHEN EVERY SECOND COUNTS

# EMERGENCY INSTRUCTION DEVICE

Imagine a life-threatening emergency has just taken place. You know to call 9-1-1, but did you know the system receives 183 million calls a year and is being inundated with non-emergency calls that are tying up operators? And once you get through to an operator, what do you do while you wait an average of 10 minutes for emergency medical services to arrive? How you respond and use the precious minutes before professional help arrives could make the difference between life and death. It's time to grab your SafetyMate.



School version, Model SM352

## SafetyMate Features

- Simple design and easy-to-follow verbal instructions can be used by anyone
- Reduces panic and confusion - takes the guesswork out by providing instructions based on yes/no questions
- Excellent for training or skills refresher
- Portable device can easily be brought to the victim's side
- Meets OSHA compliance standards for safety and first aid
- Homeland security protocols addresses Explosions and Threats such as Chemical, Radiation and Biological
- Verbal instructions given in either English or Spanish
- School version available for environments with children
- Internal memory can be updated by SafetyMate as emergency protocols and regulations change

## Technical Specifications

- **Rechargeable Batteries** SafetyMate comes with four AA rechargeable batteries and a Wall Mount Charging Cradle with AC adapter.
- **Low Battery Warning** provides both spoken and audible alerts, along with a flashing red LED, indicating that the batteries are in need of recharging, or replacement.
- **Size:** 8x11.5x3.5 inch. **Weight:** 2 lbs.
- **UL listed.** Complies to FCC standards for home or office. Rated CE compliance.
- **Each unit tested** during manufacturing.




SafetyMate protocols were developed in close cooperation with the American Red Cross and are based on the same information used by 911 operators.

"SafetyMate brings calm to those responding to a life threatening medical emergency situation and guides the rescuer to take appropriate action."

Debbie Leahy-Graves, Director,  
Health & Safety Services  
American Red Cross  
Orange County Chapter



In a medical emergency, call 9-1-1, then use SafetyMate to treat victim until professional help arrives.



1 Remove SafetyMate from Charging Cradle.

2 For ease of use, place SafetyMate close to victim.

3 Press appropriate Emergency Information Button.

4 SafetyMate will provide verbal instructions followed by a question.

5 Respond to the question by pressing Yes or No.

6 Follow verbal instructions.



Woman/Veteran Owned Small Business

## PROTECTION & DEVELOPMENT ENTERPRISES

P.O. BOX 9394 KNOXVILLE TN 37940  
PHONE 865 712 5626 FAX 865 577 1377

www.pdestore.com  
victoria@pdestore.com